

Prepping for LOINC/RELMA Class

Homework Tasks

- Download and install RELMA
 - www.LOINC.org
- Read (or skim) the RELMA manual
 - The manual downloads with the full install
- Create a matching spreadsheet



Running RELMA®



Slides downloaded and adapted with permission from: <https://www.scribd.com/doc/249865347/2014-12-03-Laboratory-LOINC-and-RELMA-Workshop>

Installing RELMA[®]



Installation Steps

- Requires Windows (or VM on a Mac)
- Make sure you have enough free disk space!
 - 2Gb is recommended
- From CD - Start – Run – <drive>:\RELMA\Setup
- Specify installation directory
 - Life will be easier if you accept the default
- Two database files installed
 - RELMA.MDB – LOINC Terms Database
 - LMOF3.MDB – Local Master Observation File
- Two sample files included
- Run from Start – (All) Programs – Regenstrief – RELMA

File Locations

- Database and Ancillary Files

- Windows XP = C:\Documents and Settings\All Users\Documents\RELMA\
- Windows Vista = C:\Users\Public\Documents\RELMA\
- Windows 7, 8 = C:\Users\Public\Public Documents\RELMA\

- Sample files

- Windows XP = C:\Documents and Settings\All Users\Documents\RELMA\Samples\
- Windows Vista = C:\Users\Public\Documents\RELMA\Samples\
- Windows 7, 8 = C:\Users\Public\Public Documents\RELMA\Samples\

Copyright Screen



The Regenstrief LOINC® Mapping Assistant (RELMA®)
Copyright Regenstrief Institute, Inc. All rights reserved.

Copyright Notice and License

The LOINC® codes, LOINC® table (regardless of format), LOINC® Release Notes, LOINC® Changes File, and LOINC® Users' Guide are copyright © 1995-2014, Regenstrief Institute, Inc. and the Logical Observation Identifiers Names and Codes (LOINC) Committee. All rights reserved.

The RELMA® program, RELMA® database and associated search index files (subject to the copyright above with respect to the LOINC® codes and LOINC® table included therein), RELMA® Community

Continue



Log in



The Regenstrief LOINC® Mapping Assistant (RELMA®)

Copyright Regenstrief Institute, Inc. All rights reserved.

Boost your LOINCing and support a noble cause
Become a LOINC Premium Member today
<http://members.loinc.org>

Username

smatney

Password

|

☒ Remember me

☐ Do not show on start up

[Register](#)


Log in

Cancel

Main Menu/Welcome Screen



Set Preferred Language



Select Linguistic Variant(s)

Regenstrief LOINC Mapping Assistant (RELMA) - Preferred Linguistic Variant

RELMA supports multiple linguistic variants which can be displayed in the search results screen. You may select your preferred linguistic variant from the list below. You may change this setting later using this dialog.

You may view all of the linguistic variants for a particular LOINC anytime using the link below.

Additional information about linguistic variants, including full text, is available at <http://loinc.org/international>.

NOTE - selecting a linguistic variant for searching may require the creation of one or more additional indexes. This process can take several minutes to complete, so please be patient. These indexes will need to be refreshed with each update to the LOINC data. Additionally, these index files will not be removed if RELMA is un-installed. You may safely delete these files at any time.

My preferred linguistic variant is...

Display	Search	Variant	Producer
<input type="checkbox"/>	<input type="checkbox"/>	Greek (GREECE)	Efstratia Kontaxi, MD, MSc, and Evripidis Stefanidis, MD, with technical su...
<input type="checkbox"/>	<input type="checkbox"/>	Italian (ITALY)	Consiglio Nazionale delle Ricerche
<input type="checkbox"/>	<input type="checkbox"/>	Italian (SWITZERLAND)	CUMUL, Switzerland
<input type="checkbox"/>	<input type="checkbox"/>	Korean (KOREA, REPUBLIC OF)	Korean Ministry for Health, Welfare, and Family Affairs
<input type="checkbox"/>	<input type="checkbox"/>	Portuguese (BRAZIL)	HL7 Brazil Institute
<input type="checkbox"/>	<input type="checkbox"/>	Russian (RUSSIAN FEDERATI...)	Yaroslavl State Medical Academy
<input type="checkbox"/>	<input type="checkbox"/>	Spanish (ARGENTINA)	Conceptum Medical Terminology Center
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Spanish (SPAIN)	the Clinical Laboratory Committee of SERVICIO EXTREMEÑO DE SALUD,...
<input type="checkbox"/>	<input type="checkbox"/>	Spanish (SWITZERLAND)	CUMUL, Switzerland

OK Cancel

Selecting a translation may require indexes to be built

LOINC Terms with Spanish Linguistic Variant

Regenstrief LOINC Mapping Assistant (RELMA) - Map Local Terms for LOINC Workshop OBR

FileToolsHIPAALab Auto MapperViewHelp

Welcome log inRegister

SearchMappingView All Working Set TermsHierarchy & Search LimitsPart SearchAnswer List Search

Local Term File

NextPreviousFirstLast

View: All

#: 1 of 190

Mapped to:Name:

OBR-4 Code:OBX-3 Code:Units:Sample Values:Limit to Default Specimen:

Click to add tag...

Accept or enter OBR name and/or OBX name

ALLERGY

Show WordsPropose TermClear InputsReset Limits

Local Term Details

Search

GridTree

Row	LOINC	Component	Property	Timing	System	Scale	Method	Rank	LongN...	ShortName	T
1	71375-0	(Alternaria...	Concentrac...	Punto...	Suero	Qn	Multidisco		ALLERGY	Inhala...	(Alternaria...
2	63199-4	(Anacardium...	Concentrac...	Punto...	Suero	Qn	Multidisco		ALLERGY	Food...	(Anacardiu...
3	63151-5	(Anacardium...	Concentrac...	Punto...	Suero	Qn	Multidisco		ALLERGY	Food...	(Anacardiu...
4	56331-2	(Anacardium...	Concentrac...	Punto...	Suero	Qn	Multidisco		ALLERGY	Food...	(Anacardiu...
5	71374-3	(Daucus...	Concentrac...	Punto...	Suero	Qn	Multidisco		ALLERGY	Food...	(Daucus...
6	71376-8	(Dermatophagoides...	Concentrac...	Punto...	Suero	Qn	Multidisco		ALLERGY	Inhala...	(Dermatop...
7	56781-8	(Pluma...	Concentrac...	Punto...	Suero	Ord	Multidisco		ALLERGY	Epider...	(Pluma...
8	11191-4	Abadejo (Polladius...	Concentrac...	Punto...	Suero	Qn			ALLERGY	Pollock...	Abadejo...
9	21470-0	Abadejo (Polladius...	Concentrac...	Punto...	Suero	Ord			ALLERGY	Pollock...	Abadejo...
10	41874-9	Abedul blanco...	Concentrac...	Punto...	Suero	Qn		1025	ALLERGY	White...	Abedul...
11	61227-5	Abedul blanco...	Fracción	Punto...	Suero	Qn			ALLERGY	White...	Abedul...
12	56422-9	Abedul blanco...	Concentrac...	Punto...	Suero	Qn			ALLERGY	White...	Abedul...
13	43841-6	Abedul blanco...	Concentrac...	Punto...	Suero	Ord			ALLERGY	White...	Abedul...

Truncated Text

Print Preview

View Details

Print Grid

Map

Same

Export

Configure Export

Configure Grid


Entry #: 1 of 190UnitsSpecimenMethodlessNo Common LimitsBattery3773 records found: 0.24s

All Spanish!

Setting User Preferences



Set User Preferences



Set User Preference Dialog

Regenstrief LOINC Mapping Assistant (RELMA) - User Preferences

General | Map Screen | Details Pages | Submissions | Local Terms Backup

Path to LOINC Database: C:\Users\Public\Documents\RELMA\RELMA.MDB Browse...

Path to Local Terms: C:\Users\Public\Documents\RELMA\RELMA.MDB Browse...

NOTE: You must restart RELMA after changing file locations.

When RELMA starts, it should

- ☒ Take me to the Welcome Screen
- ☐ Take me to the Mapping Screen
- ☐ Take me to the Simple Search Screen

☒ Do not show log in screen on startup

Apply OK Cancel

New Local Term Backup Tab

File Locations:
(See prior slide for location)

Startup Screen Preference

User Map Screen Preferences

The screenshot shows the 'User Preferences' dialog box for the Regenstrief LOINC Mapping Assistant (RELMA). The 'Map Screen' tab is selected, and several options are highlighted with callouts:


- Auto create search terms:** This option is represented by the 'Auto Select Keyword' checkbox, which is checked.
- Include/Exclude Battery Terms:** This option is represented by the 'Use Local Battery Terms in Search' checkbox, which is checked.
- Add comments on Mappings:** This option is represented by the 'Prompt for Comments When Mapping' checkbox, which is checked.

Other visible options include:


- Auto Wildcard Keywords (*):** This checkbox is unchecked.
- Display Best Search Statistics:** This checkbox is unchecked.
- Sort Search Results:** This section contains two radio buttons: 'By Component, System, Property' (selected) and 'By Score' (unchecked).
- Current Linguistic Preference:** A text field showing 'English (UNITED STATES), Regenstrief Institute Inc.' with a 'Choose...' button next to it.

At the bottom of the dialog box are three buttons: 'Apply', 'OK', and 'Cancel'.

User Detail Pages Preferences



Request Submission Preference



Regenstrief LOINC Mapping Assistant (RELMA) - User Preferences

General | Map Screen | Details Pages | Submissions | Local Terms Backup

Preferred Submission Type

☒ WEB

☐ EMAIL

Timeout in Minutes 180

Apply OK Cancel

Preferred method of request submission

New! Local Backup Preferences

Regenstrief LOINC Mapping Assistant (RELMA) - User Preferences

General | Map Screen | Details Pages | Submissions | Local Terms Backup

Automatic Backup Frequency

- ☒ Never make a backup, I'd rather do it manually.
- ☐ Make a backup every time I exit the program.
- ☐ Make a backup every days (only occurs when program is exited).

Backup Retention Policy

- ☒ Keep the last backups on disk.
- ☐ Never delete backup files. I'll clean them myself.

Apply OK Cancel

Prepare LOINC Input Spreadsheet

- In Excel – create a spreadsheet for matching with the following headers and populated with your data in the rows (sample attached):
 - Tag (optional but can assist with workflow or categories).
 - Battery Code (optional, I recommend using the form or template ID)
 - Battery Description (optional, this is the string description for the battery)
 - Local test code (required - for the observation or assessment question)
 - Description - required
 - Units of Measure (if quantitative)
- Or...you can use the attached sample spreadsheet for your mapping.



Microsoft Office
Excel Worksheet